



A Nanotech Product
The Science of Hope



High Flow Oxygen Concentrator



www.nhtech.in

MODEL: CLINOXX-NHT-10



EQUIPMENT: OXYGEN CONCENTRATOR

Model : Clintoxx-NHT-10
NET WEIGHT : 27.00KG (APPROX)
MRP : RS. 72,000/-

Distributed By:





CLINOXX-NHT-10 Oxygen concentrator device built with state of the art Electrolytic Technology using alkaline water as a primary source for generation of high purity oxygen.

CLINOXX-NHT-10 Oxygen Concentrator is a versatile cost efficient, reliable portable oxygen concentrator for your home.

Meeting National and International compliance and high-performance standards, and manufactured under Good Manufacturing Practice (GMP) implemented by AMERITTI.

CE Certified, US FDA under Good Guidance Practice and ISO 13485-2016 certified.

Every unit built with easy patient controls, bright LED that can be seen from a distance and deliver instantly (< a second) constant oxygen level (93%) and ensure patient safety and reliability for longer service intervals.

CLINOXX NHT 10 Oxygen Concentrator provides supplemental Oxygen flow therapy for patients suffering from cardiovascular disease, COPD (Chronic Obstructive Pulmonary Disease), and respiratory disease and disorders.

CLINOXX NHT 10 Oxygen Concentrator is used in residential environments, SPA/Nursing homes, Primary Health Center/Patient care facilities, etc.

Model	CLINOXX-NHT-10
Rated power (W)	160W
Energy consumption	1.4kWh/m3 of O2 under NTP
Voltage (v)	220V-240V ± 5V 50Hz ± 1Hz
Flow rate (L/min)	0.5~10
Concentration (Rated flow)	93%±5%
Sound pressure level dB (A)	≤45
Sound power level dB (A)	≤55
Net weight (kg)	27kg (Approx)
Dimension (mm)	250L×280W×420H
Alarm Signal	Alarm Signal Abnormal water level alarm; Power failure alarm; Low flow alarm.

Accessories Supplied.

Breathing Cannula.
Humidifier Bottle.
Customer Service.

NHT-10-HMT001
NHT-10-HB500
+91 7338985385



ISO-13485-2016
UK-GLOBAL-GMP
CERTIFIED

ISO-9001-2015
CE-1320 CERTIFIED